

## READ ME (2)

For subject s003 in the five-fingers task (level 3), there are two files for each repetition:

- `s003_rec_task_3_trigger_rep_Y.txt` contains the trigger information
- `s003_rec_task_3_eeg_rep_Y.txt` contains the raw EEG data

where `Y` is the repetition number.

For repetition 1, it is recommended to remove all triggers starting from trigger 101 (included) and the corresponding EEG values, based on their local timestamps.

The term `current_state` refers to the fact that participants were given the command corresponding to the current state of the hand. However, due to a software bug, the task was repeated, and the minimum number of repetitions per combination was set to 17 to ensure sufficient data.